

PRODUCT

DermaTherapy®

THE NEXT GENERATION OF HEALTHCARE TEXTILES

Prevent pressure injuries.
Improve outcomes.
Reduce costs.





DermaTherapy®

The industry's first FDA cleared, evidence-based therapeutic silk-like bedding and patient apparel to help prevent pressure ulcers.

DERMATHERAPY®:
I.V./TELEMETRY PATIENT GOWN
PILLOWCASE
FITTED SHEET
FLAT SHEET

standardtextile 

Bedding

DermaTherapy® controls the microclimate—moisture, friction and shear—between the patient’s skin and support surface. Its moisture-wicking properties keep the patient drier and cooler.

DermaTherapy® is woven with continuous filament yarns that are free of short, protruding fibers or pills that can create friction and cause irritation. The result is a silk-like sleep surface that is non-irritating and friendly to the patient’s skin.

DermaTherapy® bedding incorporates a durable antimicrobial that inhibits the survival of bacteria on the bedding. The bedding does not generate airborne particles such as lint, significantly reducing the potential for bioburden and cross-contamination.

Patient Gowns

Smooth, silk-like patient gowns move with the body, reducing friction and shear.

DermaTherapy® allows the free flow of air and moisture vapor through the fabric, enhancing patient comfort by keeping the patient drier and cooler.

Patient gowns incorporate a durable antimicrobial treatment that inhibits the survival of bacteria on the gown. Gowns do not generate airborne particles such as lint, significantly reducing bioburden.

DermaTherapy® gowns are significantly more durable, colorfast, and more stain-resistant than incumbent cotton-blend gowns. Available in standard and bariatric sizes.

Underpads

The patient-facing surface of DermaTherapy® underpads is the same material used in DermaTherapy® silk-like bedding.

Since the underpad is in direct contact with the patient, the same microclimate control of moisture, friction and shear is present to minimize skin breakdown.

Available in standard and bariatric sizes, DermaTherapy® underpads’ smooth surface also enables easier patient movement and repositioning.



DERMATHERAPY® PILLOWCASE



DERMATHERAPY® FITTED SHEET



DERMATHERAPY® FLAT SHEET



PATIENT GOWNS



I.V./TELEMETRY PATIENT GOWNS




DERMATHERAPY® UNDERPAD

Fresher



A durable antimicrobial treatment helps maintain freshness and eliminates odors on the fabric caused by germs and bacteria. The soil release finish aids with the release of stains.

DermaTherapy® fabrics do not generate airborne particles such as lint and significantly reduce bioburden on bedding and gowns.



Drier

The continuous filament yarns contain micro channels that wick moisture away from the skin 40% faster. DermaTherapy® dries 60% faster than cotton bedding.



Smoother

DermaTherapy® fabrics are engineered to reduce shear and friction on the patient's skin by eliminating the protruding short fibers found in cotton bedding.



DERMATHERAPY® MAGNIFIED



COTTON FIBERS MAGNIFIED

DermaTherapy® is a registered trademark of Precision Fabrics Group, Inc.